



HEINE 3S LED HeadLight® UNPLUGGED

Compact, light and uniquely-mobile.

The HEINE 3S LED HeadLight® UNPLUGGED is an innovative, ultra-bright, compact LED headlight featuring exceptionally-bright, coaxial illumination. The ideal headlight for general practice, ENT, dental and veterinary practitioners and surgeons.

30,000 Lux of bright, white light. With a reference colour temperature of 6,000 Kelvin and an extraordinary 30,000 Lux of illumination at a working distance of 250 mm, the 3S LED HeadLight UNPLUGGED is more than twice as bright as other standard LED and conventional headlights.

Compact and light with 100% mobility. The compact design ensures clear, unobstructed coaxial illumination. The illumination head of the 3S LED weighs only 70 grams ensuring user comfort even during long periods of use.

UNPLUGGED: No cables, no connections, no limits. With the new 3S LED HeadLight UNPLUGGED with the mPack UNPLUGGED rechargeable headband-mounted battery, you are free to work anywhere. The light weight of the battery ensures no reduction in the high level of comfort.

Versatility: The 3S LED HeadLight UNPLUGGED can be used for virtually any examination. The light spot and brightness are both adjustable and the 3S LED is compatible with all HEINE HR® and HRP® binocular loupes. A video prism optics adaptor is also available.

Many of these innovations result from our close co-operation with physicians and universities. The 3S LED HeadLight UNPLUGGED is the instrument, which physicians helped to define and develop with the HEINE R&D team.